DEFLEX PEEK

- a. [Instrucciones de uso para la técnica en inyección en aro \(ES\)](#)
- b. [Instructions for use for the ring injection technique \(EN\)](#)
- c. [Instruções de utilização para a técnica de injeção em anel \(PT\)](#)
- d. [Mode d'emploi pour la technique d'injection en anneau \(FR\)](#)
- e. [Istruzioni per l'uso per la tecnica di iniezione dell'anello \(IT\)](#)
- f. [Инструкция по применению кольцевой инъекционной техники \(RU\)](#)
- g. [Halka enjeksiyon tekniği için kullanım talimatları \(TR\)](#)
- h. [Uputstvo za upotrebu za tehniku ubrizgavanja prstena \(SR\)](#)
- i. [Instrukcja użycia dla techniki wtrysku pierścieniowego \(PL\)](#)
- j. [Návod k použití pro techniku kroužkové injekce \(CZ\)](#)
- k. [Gebrauchsanweisung für die Ringinjektionstechnik \(DE\)](#)

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Lanús (C.P.:1824) - Buenos Aires - Argentina

Tel.: (+5411) 4241-7762

info@deflex.com.ar www.deflex.com.ar

ET 80 Rev 08
04-OCT-2023

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10-DATOS DE CONTACTO

Departamento técnico: info@deflex.com.ar

B. Instructions for use for the ring injection technique (EN)

DEFLEX PEEK

1-PRODUCT DESCRIPTION

Injectable PEEK for bridge framework, crowns and partial dentures.

2-INDICATIONS / COUNTER-INDICATIONS / LIMITATIONS

- Infra-structure of bridges and crowns
- Anatomical structure of bridges and crowns
- Telescopic crowns and screw-retained crowns
- Screw-retained fixed dentures (crowns and bridges)

Note: Bridge structures only admit up to 2 intermediate elements.

3-SYMBOLS USED

	S= Small M=Medium		
	Net weight	Diameter	
	Injection temperature		
	T1= Injection or plasticized time T2=Injection holding time		
	Air pressure		
	Unique Device Identifier		

	Lot number Expiration/expiry date Manufacturing date
	Manufacturer
	Prohibited reuse Protect from the sun's rays Keep in a cool and dry place
	Do not use if packaging is damaged or defective Follow the instructions
	Authorised representative
	Medical Device

4-COMMERCIAL PRESENTATION

MEASURES AND WEIGHTS

Small / 20g
Medium / 40g

AVAILABLE COLORS

Beige
Medium pink

Primary packaging: Plastic tubes, with moisture-absorbent material

Pre-dried

Note: weight values are nominal and may differ slightly from cartridge to cartridge.

5-CARTRIDGE HANDLING

Open the plastic tube, and remove the bag containing the moisture-absorbing material. After removing the amount of bulk material you need, replace the drying bag, and close the plastic tube tightly with its lid.

6-SAFETY INFORMATION

ALLERGENIC INFORMATION

Potential Known Allergenic Substances present in the product: None

This product should not be used for patients with a suspected allergy until allergy tests have been completed to confirm that the patient is not allergic to the materials contained in the product (product sample supplied upon request for allergy test of allergy sufferers). (S21)

PEDIATRIC USE

The use of a denture base in children should be prescribed by the Pediatric Dentist. (S22)

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user is established.

7-EXPIRY / STORAGE

Store in cool dry place. (S23)

Do not expose to the rays of the sun.

Do not use if the container is damaged or defective. (S24)

Use immediately after opening. (S25)

Expiration / expiration date / Shelf life: 24 months

8-MAKING OF THE DENTURE

8.1 WORKING CAST

In cases with implants, build the working cast through an implant transfer impression. Place analogous in the impression. It is recommended to make the models in special plaster type IV (Densita).

8.2 WAXING

You can use any type of wax or castable resin.

In case of choosing castable resins coat them with a thin layer of wax to avoid damages to the investment when casting.

We recommend the use of retention beads.

In the case of bridges, the connectors between elements must not be less than 16 mm².
The use of a strengthening shoulder is recommended to ensure the forces will be supported by the infrastructure (S28) when manufacturing infrastructures to be covered with composite.

8.3 INJECTION CONDUITS (S27)

Measures: 2mm diameter x 15mm length

8.4 PLACEMENT IN THE CRUCIBLE AND RING

Bonding to the plastic crucible

Apply a thin layer of wax inside the ring to control investment expansion.

8.5 CASTING INVESTMENT

The preferred investment is phosphate based.

Percentage of EXPANSOR LIQUID / DISTILLED WATER:

70% - 30% for parallel wall bridges

60% - 40% for conical preparations

90% - 10% (or expander 100%) for multiple anchors

Manual spatulate for 90 seconds, or mechanical mixing for 60 seconds.

Pressurize at 2 kg/cm² for 15 minutes.

8.6 HEATING OF THE RING (S30)

Follow the instructions of the casting investment manufacturer.

Note: if you follow the speed heating method, we recommend not to exceed 700° C of initial preheating of the casting oven.

8.7 PREPARING THE DEFLEX PEEK CARTRIDGE

In an empty 22mm aluminum tube, fill the amount of material to be used, and trim the excess of the tube.

Calculation of amount of Deflex PEEK to use: (S32)

$[(\text{Weight of the pattern} + \text{Weight of wax conduits}) \times 1.6] + 2 \text{ grams}$

Place the cap and close the tube by folding the excess aluminum inward.

8.8 PREPARING OF THE INJECTOR


For MAD1300, MAD1400 and MD1500 injector models, remove the metal guides for muffles.

8.9 INJECTION

- Turn on the injector at least 1 hour before injecting, holding temperature at 360°C.
- Insert the brush from the cleaning kit into the oven and clean by moving the brush up and down at medium speed.

- Place (by dropping) the supplementary plunger inside the oven (Product Code PD128A – SUPPLEMENTARY PLUNGER FOR PEEK INJECTION) and then the special coin for PEEK injection (Product Code PD139 – COIN FOR PEEK INJECTION).

- Place the special Deflex high-temperature separator cartridge around it (Product Code PD140: HIGH-TEMPERATURE PEEK RELEASE AGENT), or a similar one. Apply only a thin layer with a brush, from the middle of the tube downwards (S31).

- Place the cartridge with the material in the oven of the machine.
- Remove the casting ring from the casting oven, and place it over the ring centering disc of the injector.
- Place the upper supplement for the ring and adjust the casting ring with the injector handle.
- Start the injection cycle.

DEFLEX PEEK	
	360°C + 5°C 680°F + 7°F
 T1	15' (para aros grandes) 16' (para aros pequeños)
 T2	10'
	1 Kg/cm2 14.5 PSI

-After running times T1 and T2, eject the ring, and allow to cool the casting ring to room temperature.

*Ring centering disc, leveling plunger, rings and crucible bases are part of the kit provided by Deflex.

NOTE: Injection parameters indicated for Deflex injection machine. For other injectors, consult the manufacturer or supplier of the equipment.

8.10 RING OPENING

Recover the injected part by breaking the castable investment.

8.11 RETOUCHING AND POLISHING

Retouching the flanks: burs suitable for PEEK, and roughing stones.

Polishing: with rubbers, sandpaper, goat hair brushes and Deflex gloss paste.

Sandblasting: with 50/100 µm aluminum oxide, at 2-3 bar air pressure.

Washing: in a steaming machine, or with ultrasound machine, or clean with a soapy degreaser.

NOTE: for the subsequent steps of coating with composite resin, follow the technique of your preference.

8.12 MECHANICAL AND CHEMICAL ADHESION TECHNIQUE FOR THE BONDING OF COATING MATERIALS TO PEEK.

-Sandblast with 50/100 µm aluminum oxide, at 2-3 bar air pressure

-Wash with steaming machine (with screwed work), or ultrasound.

-Let it dry well, and apply silane. Let it evaporate.

-Apply SINGLE BOND UNIVERSAL DE 3M adhesive or similar, and light cure for 25 seconds.

-Apply CERAMAGE SHOFU pre opaque paste or similar. Light cure for 1 minute.

-Apply opaque color as required and light cure for 3min.

-Continue with the conventional method of composite coating.

9-OBSERVATIONS FOR THE USER

For professional use only.

Contact your distributor or waste management local authorities to ensure proper disposal of the waste originated during the injection process. (S35)

The technical steps outlined for the processing of this product are those indicated and / or suggested by the manufacturer, and correspond to the current state of the art developed by the manufacturer.

The indications related to the injection flask refer to the Deflex brand flask.

The use of complementary techniques, materials, and equipment for processing the product is the responsibility of the user. In cases of unsuccessful results, the manufacturer limits the rights to compensation for damages exclusively to the commercial value of its products.

Do not use toothpaste or any other abrasive product since they may scratch the surface of your denture. Do not use alcohol for cleaning.

10-CONTACT DETAILS

Technical department: info@deflex.com.ar

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