DEFLEX PEEK

- A. [Instrucciones de uso para la técnica en inyección en mufla \(ES\)](#)
- B. [Instructions for use for the muffle injection technique \(EN\)](#)
- C. [Instruções de utilização para a técnica de injeção da mufla \(PT\)](#)
- D. [Mode d'emploi pour la technique d'injection de moufle \(FR\)](#)
- E. [Istruzioni per l'uso per la tecnica di iniezione della muffola \(IT\)](#)
- F. [Инструкция по применению техники впрыскивания муфеля \(RU\)](#)
- G. [Susturucu enjeksiyon tekniği için kullanım talimatları \(TR\)](#)
- H. [Uputstvo za upotrebu za tehniku ubrizgavanja prigušivanja \(SR\)](#)
- I. [Instrukcja użycia techniki wstrzykiwania mufli \(PL\)](#)
- J. [Návod k použití pro techniku vstřikování mufle \(CZ\)](#)
- K. [Gebrauchsanweisung für die Muffelinjektionstechnik \(DE\)](#)

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Lanús (C.P.:1824) - Buenos Aires - Argentina

Tel.: (+5411) 4241-7762

info@deflex.com.ar www.deflex.com.ar

ET 47 Rev 10
8-OCT-2024

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Departamento técnico: info@deflex.com.ar

B. Instructions for use for the muffle injection technique (EN)

DEFLEX PEEK

1-PRODUCT DESCRIPTION

Injectable PEEK for bridge framework, crowns and partial dentures

2-INDICATIONS / COUNTER-INDICATIONS / LIMITATIONS

- RPD's partial frameworks
- Overdentures
- Hybrid dentures
- Bridges over metal abutments
- Implant bars
- Bars
- Prostheses over attachments
- Custom abutments

Note: Bridge structures only admit up to 2 intermediate elements.

For hybrid "Toronto-type" prostheses, they can be made with a minimum of 4 abutments, evenly distributed between hemi-arches.

3-SYMBOLS USED

	S= Small M=Medium		
	Net weight	Diameter	
	Injection temperature		
	T1= Injection or plasticized time T2=Injection holding time		
	Air pressure		
	Unique Device Identifier		

	Lot number Expiration/expiry date Manufacturing date
	Manufacturer
	Prohibited reuse Protect from the sun's rays Keep in a cool and dry place
	Do not use if packaging is damaged or defective Follow the instructions
	Authorised representative
	Medical Device

4-COMMERCIAL PRESENTATION

MEASURES AND WEIGHTS

Small / 20g
Medium / 40g

AVAILABLE COLORS

Beige
Medium pink

Primary packaging: Plastic tubes, with moisture-absorbent material

Pre-dried

Note: weight values are nominal and may differ slightly from cartridge to cartridge.

5-CARTRIDGE HANDLING

Open the plastic tube, and remove the bag containing the moisture-absorbing material. After removing the amount of bulk material you need, replace the drying bag, and close the plastic tube tightly with its lid.

6-SAFETY INFORMATION

ALLERGENIC INFORMATION

Potential Known Allergenic Substances present in the product: None

This product should not be used for patients with a suspected allergy until allergy tests have been completed to confirm that the patient is not allergic to the materials contained in the product (product sample supplied upon request for allergy test of allergy sufferers). **(S1)**

PEDIATRIC USE

The use of a denture base in children should be prescribed by the Pediatric Dentist. **(S2)**

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user is established.

7-EXPIRY / STORAGE

Store in cool dry place. **(S3)**

Do not expose to the rays of the sun.

Do not use if the container is damaged or defective. **(S4)**

Use immediately after opening. **(S5)**

Expiration / expiration date / Shelf life: 24 months

8-MAKING OF THE DENTURE

NOTE: The instructions for use vary according to the application, that is, according to the type of prosthesis to be injected.

8.1 HYBRID DENTURES WITH METAL INTERFACE (Deflex UNI-BLOCK Technique) **(To cover with composite resin, conventional resin or with our M10 or ACRILATO thermo-plastics)**

8.1.1 WORKING CAST

Build the working cast through an implant transfer impression. Place analogous in the impression. It is recommended to make the models in special plaster type IV (Densita). **(S6)**

8.1.2 MOUNTING ON ARTICULATOR


Place in articulator. Semi-adjustable articulator is recommended, with clinical records.

8.1.3 DUPLICATING THE MODEL

Transfer the analogues using the splinting method (joining the abutments together, making a box in wax and filling with minimum expansion hard plaster). Cover the chimneys of the analogues with a wax thread, ensuring that the ends of the chimneys are in contact with the upper part of the duplicating flask. Duplicate with gelatine for plaster or silicone.

After the duplication material has hardened, it is recommended to first find the wax chimneys, remove the screws, in order to then be able to separate the working cast.

Place new analogues in the duplicate and screw them in. Empty with special type IV plaster.

8.1.4 MOUNTING THE DUPLICATE MODEL

With the bite registration or taking a new registration from the already prepared working cast, mount the duplicate cast in the articulator with a new base. It is in this duplicate model where we will do the alignment and the injection.

8.1.5 PREPARATION OF ABUTMENTS

Place the metal abutments in the duplicated model. Make retentions on them, with a metal cutting disc and sandblasting with 50/100 µm aluminium oxide, at 2-3 bar air pressure.


8.1.6 ALIGNMENTS WITH ARTIFICIAL TEETH

Make the wax pattern and align the artificial teeth after making mechanical retentions (central drilling of 2.3 mm, and by proximal 1.5 mm).

8.1.7 PLACEMENT IN THE MUFFLE

Cover the chimneys of the analogues with putty or hard silicone to the edge of the abutment. Use special type IV plaster.

8.1.8 INJECTION CONDUITS (S11)

Perform tubes in wax. The tubes must follow a straight path on the material income plane during the injection. For this reason, the model must be placed on that same plane; it must neither be buried nor be much higher than the plane of the muffle. The conduits are positioned and glued on the abutments. Diameter of the ducts: 5-6 mm, in the shape of a half round.

8.1.9 COUNTER MUFFLE

Cover the chimneys of the remaining modelling with the same plaster of the counter muffle. Apply plaster in two stages, using plaster separator after first coat of plaster has set. Use special type IV plaster.

8.1.10 MUFFLE CLEANSING

Loosen the screws two turns before cleansing. Immersion time: from 6 to 8 minutes. Degrease with soapy water.

8.1.11 TUBES CORRECTION

Rectify and correct all the tubes, eliminating any thin layer or excess that may have remained. *Note: for the peek injection, no type of separator is placed.*

8.1.12 INJECTION PARAMETERS

Turn on the injector 30 minutes before the injection with the following parameters.

DEFLEX PEEK	
	390°C/395°C
 T1	15'
 T2	2'
	6 Kg/cm2 +1 86 PSI + 14

NOTE: Injection parameters indicated for Deflex injection machine. For other injectors, consult the manufacturer or supplier of the equipment.

8.1.13 HEATING THE MUFFLE

Immerse the closed muffle in boiling water for 20 minutes.

8.1.14 PREPARING THE DEFLEX PEEK CARTRIDGE

In an empty 22mm aluminium tube, fill the amount of material to be used. Cut the excess aluminium from the tube.

Superior hybrids dentures: 28 grams

Lower hybrids dentures: 30 grams

NOTE: The amounts indicated are only suggested values based on user experience. It can vary significantly depending on the dimensions of the prosthesis to be injected. Using less material than what the prosthesis to inject requires will imply the repetition of all the work.

Another option is to weigh all the wax used (pattern of the modelling and ducts) and multiply by 1.6, and to that result add 3 gr.

Place the cap and close the tube by folding the excess aluminum inward.


8.1.15 INJECTION

Insert the brush from the cleaning kit into the oven and clean by moving the brush up and down at medium speed.


Place (by dropping) the supplementary plunger inside the oven (Product Code PD128A – SUPPLEMENTARY PLUNGER FOR PEEK INJECTION) and then the special coin for PEEK injection (Product Code PD139 – COIN FOR PEEK INJECTION).

Place the special Deflex high-temperature separator cartridge around it (Product Code PD140: HIGH-TEMPERATURE PEEK RELEASE AGENT), or a similar product. Apply only a thin layer with a brush, from the middle of the tube downwards.

Place the tube with DEFLEX PEEK in the injection oven, position the mold - after draining it well upon removing it from the water in which it was preheated - and start the injection cycle.

8.1.16 MUFFLE OPENING

Eject the muffle from the injector and allow to cool it down.

NOTE: Once the injection cycle is finished, it is very important to clean the injector oven with the Deflex cleaning brush and a drilling or screwing machine. Product code: PD138 (PEEK CLEANING BRUSH)

For the cleaning process, the injector must be at 385°C.

De-flask by gently tapping the plaster on top of the counter-flask.

8.1.17 RETOUCHING AND POLISHING

Retouching the flanks: burs suitable for PEEK, and roughing stones.

Polishing: with rubbers, sandpaper, goat hairbrushes and Deflex gloss paste.

Sandblasting: with 50/100 µm aluminium oxide, at 2-3 bar air pressure.

Washing: in a steaming machine with the screws on, or with ultrasound machine, or clean with a soapy degreaser.

NOTE: for the subsequent steps of coating with composite resin, follow the technique of your preference.

8.2 BRIDGES

Infra-structures to be covered, or structures with an anatomical finish.

8.2.1 PREPARATION OF THE WORKING MODEL AND ASSEMBLY ON THE ARTICULATOR.

We will build a working model and mount it conventionally on an articulator.

8.2.2 METALLIC ABUTMENTS

Place the metallic abutments on the working model. Make mechanical retentions on them, and sandblast with 50 µm aluminium oxide, at 3-4 bar air pressure.

8.2.3 WAXING

Make the wax pattern. The connectors between elements must not be less than 16 mm².

Note: Now you can choose to build the pattern, lift it and flask it, or do the technique described above (duplicating the model) and its further development (see 8.1.3.).

Note: with this technique of model duplication and injection on metallic abutments it is also possible to make bars of different types, and to customize abutments for metal-free crown work.

8.2.4 PLACEMENT IN THE MUFFLE

Lift the pattern from the model, screw on analogues, and place in the flask with special type IV plaster. Bury up to the third section of the cervical.

8.2.5 INJECTION CONDUITS (S11)

Perform tubes in wax. The wax conduits must follow a straight path on the material income plane during the injection. For this reason, the model must be placed on that same plane; it must neither be buried nor be much higher than the muffle plane.

Cover the chimneys of the analogues with hard putty. Apply plaster separator.

Diameter of the ducts: 4 mm.

8.2.6 COUNTER MUFFLE

Cover the chimneys of the abutments with hard silicone mass.

Apply plaster in two stages, using plaster separator after first coat of plaster has set. Use special type IV plaster.

8.2.7 MUFFLE CLEANSING

Remove the screws and submerge in boiling water. Immersion time: from 6 to 8 minutes.

Degrease with soapy water.

8.2.8 TUBES CORRECTION

Rectify and correct all the tubes, eliminating any thin layer or excess that may have remained.

Note: for the peek injection, no type of separator is placed.

8.2.9 INJECTOR PREPARATION - INJECTION PARAMETERS - MUFFLE HEATING - CARTRIDGE PREPARATION - INJECTION - MUFFLE OPENING - TOUCHING AND POLISHING

These technical steps are the same as those already seen from point 8.1.12 to 8.1.17, except for the amount of PEEK to use, which must follow this formula:

[(Weight of the pattern + Weight of wax conduits) x 1.6] + 3 grams

8.2.10 MECHANICAL AND CHEMICAL ADHESION TECHNIQUE FOR THE BONDING OF COATING MATERIALS TO PEEK.

- Sandblast with 50/100 µm aluminium oxide, at 2-3 bar air pressure
- Wash with steaming machine (with screwed work), or ultrasound.
- Let it dry well and apply silane. Let it evaporate.
- Apply SINGLE BOND UNIVERSAL DE 3M adhesive, and light cure for 25 seconds.
- Apply CERAMAGE SHOFU pre opaque paste or similar. Light cure for 1 minute.
- Apply opaque colour as required and light cure for 3min.
- Continue with the conventional method of composite coating.

8.3 RPD's FRAMEWORKS

(To coat with other injectable material or conventional acrylic resin)

8.3.1 MODEL CASTING

It is recommended to make the models in special plaster type IV.

8.3.2 PARALLELIZATION OF THE MODEL (S7)

In the same way as with chrome, it is recommended to use a pendular parallelometer to look for the most favourable insertion axle, determining the dental equator, according to which the retainers we will be designed.

Then, proceed to eliminate all the inconvenient or blind spots, applying wax and passing the rod.

8.3.3 RETENTIVITY (S8)

Once the parallelization has been completed, tidy up any wax excess that may have been left over the dental equator and the gingival pad, from the medial area until the retainer's end. This will give the dentures the necessary retentivity.

8.3.4 RELIEFS

Perform reliefs in wax covering the cervical areas.

8.3.5 DUPLICATE

Make a duplicate of the model to work on hereafter.

IMPORTANT! Perform the casting with plaster type IV (Densita), which is able to withstand the pressure to be used in the injection.

8.3.6 DEFINITIVE WAXING

A wax pattern with the desired shape is made on the duplicated model.

Palatal plate: 1 to 1.5 mm

Lingual bar: not less than 2.5 mm thick and as wide as possible in the shape of a half pear (wide part towards the floor of the mouth).

Opposing retainer: 1.5mm thick in the rigid third tapering towards the end and 3mm wide in the shape of a half round.

Active retainer: rigid third 1.3 / 1.5 mm thick tapering towards the end to at least 1 mm, and in the width of 2.5 mm from the rigid third sector tapering towards the elastic third.

Occlusal rests: 0.8 / 1 mm minimum, and its connector not less than 1 mm.

Mesh: thickness not less than 1 mm thick.


8.3.7 PLACING IN THE MUFFLE

Scrub the retentive areas of the plaster teeth.

Use special plaster type IV

8.3.8 INJECTION CONDUITS (S11)

Perform tubes in wax of 6mm diameter. The tubes must follow a straight path on the material income plane during the injection. For this reason, the model must be placed on that same plane; it must neither be buried nor be much higher than the muffle plane. Direct the lateral main sprues towards molar / premolar, in a "V" shape, and not in a "T" shape.

8.3.9 COUNTER MUFFLE

Put plaster separator. Present counter muffle and screw. Fill with 250 g of special type IV plaster, in two stages (100 g and 150 g, respectively). Let harden for at least 30 minutes.

8.3.10 MUFFLE CLEANSING

Remove the screws and submerge in boiling water. Immersion time: from 6 to 8 minutes. Degrease with soapy water.

8.3.11 TUBES CORRECTION

Rectify and correct all the tubes, eliminating any thin layer or excess that may have remained.

8.3.12 INJECTION PARAMETERS

Turn on the injector 30 minutes before the injection with the following parameters.

DEFLEX PEEK	
	390°C/395°C
 T1	15'
 T2	2'
	6 Kg/cm ² +1 86 PSI + 14

NOTE: Injection parameters indicated for Deflex injection machine. For other injectors, consult the manufacturer or supplier of the equipment.

8.3.13 HEATING THE MUFFLE

Immerse the closed muffle in boiling water for 20 minutes.

8.3.14 PREPARING THE DEFLEX PEEK CARTRIDGE

In an empty 22mm aluminium tube, fill the amount of material to be used. Cut the excess aluminium from the tube.

Calculation of amount of Deflex PEEK to use:

[(Weight of the pattern + Weight of wax conduits) x 1.6] + 3 grams

Place the cap and close the tube by folding the excess aluminum inward.


8.3.15 INJECTION

Insert the brush from the cleaning kit into the oven and clean by moving the brush up and down at medium speed.


Place (by dropping) the supplementary plunger inside the oven (Product Code PD128A – SUPPLEMENTARY PLUNGER FOR PEEK INJECTION) and then the special coin for PEEK injection (Product Code PD139 – COIN FOR PEEK INJECTION).

Place the special Deflex high-temperature separator cartridge around it (Product Code PD140: HIGH-TEMPERATURE PEEK RELEASE AGENT), or a similar product. Apply only a thin layer with a brush, from the middle of the tube downwards.

Place the tube with DEFLEX PEEK in the injection oven, position the mold - after draining it well upon removing it from the water in which it was preheated - and start the injection cycle.

8.3.16 MUFFLE OPENING

Eject the muffle from the injector and allow to cool it down.

NOTE: Once the injection cycle is finished, it is very important to clean the injector oven with the Deflex cleaning brush and a drilling or screwing machine. Product code: PD138 (PEEK CLEANING BRUSH)

For the cleaning process, the injector must be at 385°C.

De-flask by gently tapping the plaster on top of the counter-flask.

8.3.17 RETOUCHING AND POLISHING

Retouching the flanks: burs suitable for PEEK, and roughing stones.

Polishing: with rubbers, sandpaper, goat hairbrushes and Deflex gloss paste.

Washing: in a steaming machine with the screws on, or with ultrasound machine, or clean with a soapy degreaser.

NOTE: for the subsequent steps of coating with composite resin, follow the technique of your preference.

NOTE: the teeth are assembled according to the Deflex injectable material technique, or it can also be done with the conventional technique for acrylic.

(S9)

9-OBSERVATIONS FOR THE USER

NOTE: We recommend cleaning the oven with the DEFLEX CLEANING KIT. For more information, consult the instructions for use of the CLEANING KIT (ET91).

For professional use only.

Contact your distributor or waste management local authorities to ensure proper disposal of the waste originated during the injection process. **(S35)**

The technical steps outlined for the processing of this product are those indicated and / or suggested by the manufacturer, and correspond to the current state of the art developed by the manufacturer.

The indications related to the injection flask refer to the Deflex brand flask.

The use of complementary techniques, materials, and equipment for processing the product is the responsibility of the user. In cases of unsuccessful results, the manufacturer limits the rights to compensation for damages exclusively to the commercial value of its products.

Do not use toothpaste or any other abrasive product since they may scratch the surface of your denture. Do not use alcohol for cleaning.

10- CONTACT DETAILS

Technical department: info@deflex.com.ar

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